Overview 0000 Measuring Single Genes

Measuring Multiple Genes

Lecture 2: Early Transcriptomic Strategies BIOINF3005/7160: Transcriptomics Applications

Dr Stephen Pederson

Bioinformatics Hub, The University of Adelaide

March 16th, 2020



Overview 0000 Measuring Single Genes

Measuring Multiple Genes

Overview

Measuring Single Genes

Measuring Multiple Genes



Overview •000 Measuring Single Genes

Measuring Multiple Genes

Overview



Measuring Multiple Genes

The Motivation

- The transcriptome is a highly dynamic set of molecules
- Small changes can potentially have significant ramifications
 - e.g. a "Master Regulator" can determine cellular fate
- RNA molecules are small
 - How do we find what's in our sample?
 - How do we quantify RNA?
 - And how do we compare one or more groups?



Overview 0000 Measuring Single Genes

Measuring Multiple Genes

Technological Developments

- Technological developments are constant
- Technologies are often transient
- Key technologies are:
 - 1. Real Time Polymerase Chain Reaction (RT-PCR)
 - 2. Expressed Sequence Tags (EST)
 - 3. Serial/Cap Analysis of Gene Expression (SAGE/CAGE)
 - 4. Microarray technologies
 - 5. Sequencing technologies
- Analytic methodologies often lag technologies



Measuring Multiple Genes

A Simplified History



EST (blue); SAGE / CAGE (yellow); Microarrays (red); RNA Seq (black)¹



¹Rohan Lowe et al. "Transcriptomics technologies". In: *PLOS Computational Biology* 13.5 (May 2017), pp. 1–23. DOI: 10.1371/journal.pcbi.1005457. URL: https://doi.org/10.1371/journal.pcbi.1005457.

Overview 0000 Measuring Single Genes

Measuring Multiple Genes

Measuring Single Genes



Measuring Multiple Genes

The Northern Blot

- One of the earliest strategies²
- Developed as an extension of the Southern Blot³ (DNA)
- Gel Electrophoresis-based strategy
 - · Based on size differentiation and probe sequences

³E. M. Southern. "Detection of specific sequences among DNA fragments separated by gel electrophoresis". In: J. Mol. Biol. 98.3 (1975), THE UNIVERSITY #ADELAIDE

²J. C. Alwine, D. J. Kemp, and G. R. Stark. "Method for detection of specific RNAs in agarose gels by transfer to diazobenzyloxymethyl-paper and hybridization with DNA probes". In: *Proc. Natl. Acad. Sci. U.S.A.* 74.12 (1977), pp. 5350–5354.

Measuring Multiple Genes

The Northern Blot

- RNA is extracted then denatured
- RNA is size separated using Gel Electrophoresis
- RNA is transferred to a "blotting membrane"
- Treat the membrane with a labelled probe
 - Probes are complementary to the "target sequence"
 - Probes are labelled with fluorescent dye or radioactive atoms



Measuring Multiple Genes

The Northern Blot



Taken from Ramkumar et al., "Effect of orientation of transcription of a gene in an inverted transferred DNA repeat on transcriptional geneTHE UNIVERSITY silencing in rice transgenics-a case study"

Measuring Multiple Genes

The Northern Blot

- Prominent usage before genomes were sequenced
- Can possibly detect different isoforms
- Crude quantitation using Densitometric Analysis
 - What limitations might this have?



Measuring Multiple Genes

RT-qPCR

- Reverse Transcriptase quantitative PCR
 - Sometimes called: qPCR, RT-PCR
- Often considered to be the "gold standard" for quantitation
- Targets a specific transcribed region via specific primers
 - Primers must be individually designed
 - Primers often span exon-exon junctions



Measuring Multiple Genes

RT-qPCR

- 1. Reverse Transcriptase converts RNA to cDNA
 - Primers are required: Can target poly-A or random
- 2. Sequence-specific primers amplify the target fragment in cycles
 - Fluorescent dye is commonly incorporated during amplification
- 3. Abundance of target will grow exponentially (\times 2) for each amplification cycle
- 4. The cycle where abundance reaches the "limit of detection" is estimated (C_T)



Measuring Multiple Genes

RT-qPCR





Modified from By Lokeshthimmana - Own work; CC BY-SA 4.0; https://commons.wikimedia.org/w/index.php?curid=76313637

Overview 0000 Measuring Single Genes

Measuring Multiple Genes





©Lokesh Thimmana, under the guidance of Dr. G. Mallikarjuna, Assistant Professor, Molecular Biology, Agri Biotech Foundation.



Modified from By Lokeshthimmana - Own work; CC BY-SA 4.0; https://commons.wikimedia.org/w/index.php?curid=76313637

RT-qPCR





⁴Ma Mingxiao et al. "TaqMan MGB Probe Fluorescence Real-Time Quantitative PCR for Rapid Detection of Chinese Sacbrood Virus". In: MADELAIDE

Measuring Multiple Genes

```
RT-qPCR
```

- Can be used with a standard curve and dilution series to estimate absolute quantity of an RNA *within a sample*
- Can be used to compare across samples for relative abundance



Measuring Multiple Genes

```
RT-qPCR
```

- Can be used with a standard curve and dilution series to estimate absolute quantity of an RNA *within a sample*
- Can be used to compare across samples for relative abundance

What may be a fundamental issue when comparing across samples?



Measuring Multiple Genes

Normalisation

- There may be pipetting and other technical differences between samples
 - These are **non-biological** in origin
- To correct for these we can normalise our data
- In RT-qPCR this is often done using "housekeeper" genes
 - We choose genes which *should not* change between samples/groups
 - These are commonly structural genes such as $ACTN\beta$ or GAPDH



Measuring Multiple Genes

Estimating Change In Expression

- Relative abundances are often referred to as fold-change (FC)
 - Down regulation is squeezed between 0 and 1
 - Up regulation ranges from 1 to ∞
- We often use log₂ fold-change to get a better scale, e.g.
 - A 2-fold increase in abundance: $\log_2 2^1 = 1$
 - A 2-fold decrease in abundance: $\log_2 \frac{1}{2} = \log_2 2^{-1} = -1$
 - No change in abundance $\log_2 1 = \log_2 2^0 = 0$
- This is often abbreviated as *logFC*



Measuring Multiple Genes

Estimating Change In Expression

- For RT-qPCR the estimate of logFC is known as $\Delta\Delta C_T$
- To calculate this, we calculate **two** changes in C_T
 - 1. ΔC_T relative to the housekeeper(s)
 - 2. $\Delta\Delta C_T$ across samples for our gene/fragment of interest
- The first step corrects for technical errors
- The second step estimates our true change in abundance



Overview 0000 Measuring Single Genes

Measuring Multiple Genes

Estimating Change In Expression

Within each sample

$$\Delta C_{T} = C_{t[\text{gene}]} - C_{t[\text{HK}]}$$

Across samples/groups

$$\Delta \Delta C_{\mathcal{T}} = -(\Delta C_{\mathcal{T}[\text{group1}]} - \Delta C_{\mathcal{T}[\text{group2}]})$$

This formulation assumes equal amplification efficiency for all primers/genes (i.e. Efficiency = 2)



Measuring Multiple Genes

Estimating Change In Expression

- Housekeeper genes must be matched to the "gene of interest" within each sample and within each qPCR reaction
- Choosing > 1 housekeeper gene is advised
- Measurements are often taken in triplicate/quadruplicate for each sample (reactions sometimes fail)

Both Northern blots and RT-qPCR use targeted primers, but in very different ways



Measuring Multiple Genes

Measuring Multiple Genes



Measuring Multiple Genes

Expressed Sequence Tags

- The first attempt at capturing the larger transcriptome was via Expressed Sequence Tags⁵ (ESTs) in 1991
 - Sequenced 609 mRNA human brain mRNA sequences
 - ESTs were generated by reverse transcribing poly-A selected mRNA, amplified using random primers
 - Used ESTs $\sim 100-800$ nt
 - Obtained actual sequences using Sanger Sequencing
- >10 years before the Human Genome Project completed
- Just discovering genes was a huge priority



⁵Mark D. Adams et al. "Complementary DNA Sequencing: Expressed Sequence Tags and Human Genome Project". In: Science 252.5013 (1995 UNIVERSITY pp. 1651–1656. ISSN: 00368075, 10959203. URL: http://www.jstor.org/stable/2876333.

Measuring Multiple Genes

Expressed Sequence Tags





Source: https://www.khanacademy.org/science/high-school-biology/hs-molecular-genetics/hs-biotechnology/a/dna-sequencing

Measuring Multiple Genes

Serial Analysis of Gene Expression

Serial Analysis of Gene Expression⁶ (SAGE) was the first attempt to quantify expression on a larger scale

- 1. Conversion of mRNA to ds-cDNA using biotinylated primers (often poly-T)
- 2. cDNA is bound to beads using biotin and cleaved
- 3. 11-mer "tags" were produced after cleavage and concatenated
- 4. Sequenced by Sanger Sequencing
- 5. Tags were "de-convoluted" and counted



⁶V. E. Velculescu et al. "Serial analysis of gene expression". In: Science 270.5235 (1995), pp. 484-487.

Measuring Multiple Genes

Serial Analysis of Gene Expression





Also see: https://www.scq.ubc.ca/wp-content/uploads/2006/07/SAGE3b.gif

Measuring Multiple Genes

Serial Analysis of Gene Expression

- The word "tag" is still commonly used in some NGS manuals and software
- The term "Digital Gene Expression" arose during this era
 - Is sometimes shortened to DGE, but **does not** stand for *Differential* Gene Expression.
- SAGE doesn't rely on probes targeting known sequences
- Variants on the technique are still used⁷
 - Even used these concatenated tags in early NGS contexts⁸

⁷A. M. Zawada et al. "Massive analysis of cDNA Ends (MACE) and miRNA expression profiling identifies proatherogenic pathways in chronic kidney disease". In: *Epigenetics* 9.1 (2014), pp. 161–172.



⁸H. Matsumura et al. "SuperSAGE array: the direct use of 26-base-pair transcript tags in oligonucleotide arrays". In: Nat. Methods 3.6 (2000) UNIVERSITY #ADELAIDE

Measuring Multiple Genes

Cap Analysis of Gene Expression

- A variant technique is Cap Analysis of Gene Expression⁹
- Targets Transcription Start Site (TSS) of mRNA via the 5' cap
 - Specifically for identification of the exact TSS and analysis of promoters
- Original 27nt long, but now only limited by NGS length
- Heavily used in FANTOM (Functional ANnoTation Of the Mammalian genome) project



⁹R. Kodzius et al. "CAGE: cap analysis of gene expression". In: Nat. Methods 3.3 (2006), pp. 211–222.

Measuring Multiple Genes

SAGE Vs CAGE

• Primers which target the poly-A sequence will capture mature mRNA

- mRNA will also be intact (i.e. not degraded)
- CAGE targets transcriptional initiation
 - Transcripts may not be "mature"
 - 5' Cap must be in place (i.e. not degraded)
- Both techniques still involve concatenation of "tags"



Measuring Multiple Genes

Microarrays



EST (blue); SAGE / CAGE (yellow); Microarrays (red); RNA Seq (black)¹⁰



¹⁰Rohan Lowe et al. "Transcriptomics technologies". In: PLOS Computational Biology 13.5 (May 2017), pp. 1–23. DOI: 10.1371/journal.pcbi.1005457. URL: https://doi.org/10.1371/journal.pcbi.1005457.